

SleepAngel®

Improved patient and caregiver safety with SleepAngel Barrier bedding – barrier for viruses, allergens, pathogens and liquids

Did you know, The pillow (and soft surfaces overall) is an oft-overlooked link in a chain of cross-infection

Traditional pillows require add-on pillow protectors & pillow covers. This increases your purchase cost & ongoing laundry costs

Introducing SleepAngel Barrier Bedding for the Healthcare industry to address the post COVID-19 challenges and fight Healthcare Acquired Infections

SleepAngel PneumaPure™ technology, is globally patented (incl. EU, UK, North America, Japan, China, Russia, India, Middle East). Products are certified as Class I Medical Devices in the EU and registered Class I Medical Device with the FDA.



Durable, soft touch and high performance textile that blocks out contaminants.* Vapour permeable and waterproof.

Heat sealed seams provide complete barrier - no stitching holes.

Breathable mechanical filter that allows clean air flow into the pillow but blocks out pathogens and allergens.



Certified as Class I Medical devices in the EU
Registered Class I Medical Device with the FDA

* Tested pathogens include
Coronavirus 229E (HcoV 229E), MRSA, E.Coli, Adenovirus 5, Pseudomonas aeruginosa, C.Difficile, Candida albicans, Aspergillus niger, Norovirus, Influenza Type A

	SleepAngel Pillow	Traditional Pillow
User comfort	Yes	Yes
Protected against Human Coronavirus	Yes	No
Asthma and Allergy friendly	Yes	N/A
Always stays clean inside with barrier for Viruses*, Allergens, Bacteria & Liquids	Yes	N/A
Ease of cleaning & maintenance	Yes, just with wipes	N/A
Periodic laundry wash	Not Required	Required
Warranty	2 Years	N/A

SleepAngel®



SleepAngel Duvet

Fabric face: 100% Polyurethane, Royal blue or Medical white
Backing: 100% Polyamide
Fabric weight: 180 g/m²
Filling: 100% polyester thermobonded wadding 200g/m²
Label: Print (White print on royal blue fabric;
light grey print on medical white fabric)



SleepAngel Mattress

Fabric face: 100% Polyurethane, Dark blue and Antislip gray
Backing: 100% Polyamide
Fabric weight: 230 g/m²
Bottom: Anti slip
Label: Print (White print on dark blue fabric)
Filling: Low, medium-high and high risk.
Top: Castellated soft memory foam. Gel inserted on heel area.
Middle part: Soft high resilience polyurethane foam.
Side walls and base: Hard high resilience polyurethane foam.



SleepAngel Positioners

Shapes: SOAP, TRIANGULAR, WEDGE, CYLINDER
Outer cover: Outer fabric is a multi-stretch polyurethane coated textile. Vapour permeable and waterproof.
Easy to wipe clean using approved disinfectants.
Fabric face: 100% Polyurethane, Royal blue or Medical white
Backing: 100% Polyamide
Fabric weight: 180 g/m²
Filling: Different foam types (e.g. memory foam)
Label: Print (White print on royal blue fabric;
light grey print on medical white fabric)



Clinically proven barrier to pathogen ingress, yet well ventilated, breathable and comfortable.



Patented bedding range, clinically proven to block viruses, allergens, microbes and mold.



Hermetically sealed cushioning device keeps pathogens out, because there is no stitching.